



TECHNICAL DATA

Components	up to 96
Valves' Diameter	DN06
Dispensing Accuracy	Up to 0.01 g
Dispensing Time	4 min for a 2.5 kg formula using 4 components



Campione

Laboratory Units

ADVANTAGES

Extremely easy to operate with user friendly interface which can be linked with the most common ERP and MRP programmes

Easy interface with any brand of spectrophotometer

Consistent formula reproduction insuring perfect product repeatability

Single point dispensing; it can dispense into containers of any size

Elimination of wastage

Complete tracking of operation and ISO compliance traceability

Password protected formula database and automatic back-up to prevent data loss

Production planning

Easy maintenance

Improved workplace cleanliness, tidiness and safety

This system automatically produces laboratory samples for QC with the highest degree of accuracy and reproducibility

Thanks to its **center point dispensing technology**, the CAMPIONE can dispense into containers of any size

The CAMPIONE is normally combined with a larger production system sharing the same components in order to **ensure absolute consistency between QC and production batches**

The CAMPIONE can also be **installed as a stand-alone system** having its own pumping structure whenever the

production of highly accurate small batches is required.

Wet washing cycle **starts automatically at the end of every production.**

The size of the dispensing head and the number of valves **can be adapted following customer's requirements.**

The Campione can handle both **water based as well as solvent-based products** (Ex-ATEX or IEC-EX - version)